

HEATING BLOCK

LCD display. It easy to setup and use

Simultaneous display of set and actual time, temperature

Temperature deflection adjusting meet with different demands

It has to equipped with various optional blocks, and supplied with these blocks mentioned below.

Block for 0.2ml Tubes

Block for 0.5ml Tubes

Block for 1.5 and 2.0mL tubes

The instrument can adapt to different tubes & wells to cope with experiments needs.

Easily replicable metal blocks and is very simple to clean and sterilize

Microprocessor controlled incubation time and temperature

Beep-signal / Stop after program completion

Peltier design of MK-20 provides efficient thermal control

Maximum temperature 100 degree Celsius

Minimum temperature is ambient temperature.