

- Superior quality—stringent quality control and industry leading manufacturing process
- Convenient color-coded Five Buffer System
- Includes universal Tango buffer for double-digestions
- BSA premixed in reaction buffers
- Wide selection of restriction endonuclease specificities
- 100% activity of all FastDigest enzymes in the universal buffer
- 100% buffer compatibility with downstream applications
- Complete digestion in 5–15 minutes
- Direct loading on gels
- No star activity

Quantity: 250µL

Product Line: FastDigest

Methylation Sensitivity: Not dam methylation-sensitive, Not dcm methylation-sensitive, CpG methylation-sensitive

Enzyme: Not I

Compatible Buffer: 10x FastDigest Buffer/FastDigest Green Buffer

Sensitive to Heat Inactivation: Yes

Optimal Reaction Temperature: 37° C

Type IIS RE: No